



Vector-borne
Infections

AmpliRun[®] PCR Controls

The best tool for NAT processes



Obtain reliable &
comparable results




Validate and verify
new assays



Comply with lab
accreditation requirements

- Independent external run controls.
- Purified complete microbial genome.
- Any sequence can be amplified.
- Concentration range determined by qPCR.
- Suitable for all molecular testing platforms.
- Non infectious.
- Lyophilized presentation.



vircell 
MICROBIOLOGISTS

AmpliRun® Amplification Controls

Kit contents:

- 1 vial with highly concentrated DNA/RNA.
- 1 vial with molecular biology grade water.

Description	Cat. No.
AMPLIRUN® BARTONELLA QUINTANA DNA CONTROL	MBC006
AMPLIRUN® BORRELIA AFZELII DNA CONTROL	MBC078
AMPLIRUN® BORRELIA BURGENDORFERI DNA CONTROL	MBC076
AMPLIRUN® BORRELIA GARINII DNA CONTROL	MBC077
AMPLIRUN® CHIKUNGUNYA VIRUS RNA CONTROL	MBC099
AMPLIRUN® DENGUE 1 VIRUS RNA CONTROL	MBC055
AMPLIRUN® DENGUE 2 VIRUS RNA CONTROL	MBC056
AMPLIRUN® DENGUE 3 VIRUS RNA CONTROL	MBC057
AMPLIRUN® DENGUE 4 VIRUS RNA CONTROL	MBC058
AMPLIRUN® EASTERN EQUINE ENCEPHALITIS RNA CONTROL	MBC097
AMPLIRUN® FRANCISELLA TULARENSIS DNA CONTROL	MBC110
AMPLIRUN® LEISHMANIA CHAGASI DNA CONTROL	MBC060
AMPLIRUN® LEISHMANIA INFANTUM DNA CONTROL	MBC032
AMPLIRUN® ORIENTIA TSUTSUGAMUSHI DNA CONTROL	MBC106
AMPLIRUN® RICKETTSIA CONORII DNA CONTROL	MBC042
AMPLIRUN® ROSS RIVER VIRUS RNA CONTROL	MBC130
AMPLIRUN® ST LOUIS ENCEPHALITIS VIRUS RNA CONTROL	MBC101
AMPLIRUN® TICK BORNE ENCEPHALITIS VIRUS RNA CONTROL	MBC045
AMPLIRUN® TOSCANA VIRUS RNA CONTROL	MBC129
AMPLIRUN® TRYPANOSOMA CRUZI DNA CONTROL	MBC059
AMPLIRUN® TRYPANOSOMA RANGELI DNA CONTROL	MBC093
AMPLIRUN® VENEZUELAN EQUINE ENCEPHALITIS RNA CONTROL	MBC096
AMPLIRUN® WEST NILE VIRUS RNA CONTROL	MBC069
AMPLIRUN® WESTERN EQUINE ENCEPHALITIS RNA CONTROL	MBC098
AMPLIRUN® YELLOW FEVER VIRUS RNA CONTROL	MBC100
AMPLIRUN® ZIKA VIRUS (ASIAN LINEAGE) RNA CONTROL	MBC103
AMPLIRUN® ZIKA VIRUS RNA CONTROL	MBC072
OTHER	
AMPLIRUN® BARTONELLA HENSELAE DNA CONTROL	MBC005
AMPLIRUN® BRUCELLA ABORTUS DNA CONTROL	MBC009
AMPLIRUN® COXIELLA BURNETII DNA CONTROL	MBC018
AMPLIRUN® TOXOPLASMA GONDII DNA CONTROL	MBC047

AmpliRun® Total Extraction & Amplification Controls

Kit contents:

- 10 vials for single use that contain the complete microorganism in a matrix that mimics human specimens.

Description	Cat. No.
AMPLIRUN® TOTAL ZIKV-DENV -CHIKV CONTROL (PLASMA) Zika virus, dengue virus (1-4), chikungunya virus	MBTC023