



Sexually
Transmitted
Infections

AmpliRun® PCR Controls

The best tool for NAT processes



Obtain reliable & comparable results




Validate and verify new assays



Comply with lab accreditation requirements

- Independent external run controls.
- Purified complete microbial genome.
- Any sequence can be amplified.
- Concentration range determined by qPCR.
- Suitable for all molecular testing platforms.
- Non infectious.
- Lyophilized presentation.



vircell 
MICROBIOLOGISTS

AmpliRun® Amplification Controls

RUO

Kit contents:

1 vial with highly concentrated DNA/RNA.

1 vial with molecular biology grade water.

Description	Cat. No.
AMPLIRUN® CHLAMYDIA TRACHOMATIS DNA CONTROL	MBC012
AMPLIRUN® HAEMOPHILUS DUCREYI DNA CONTROL	MBC021
AMPLIRUN® HERPES SIMPLEX 1 DNA CONTROL	MBC023
AMPLIRUN® HERPES SIMPLEX 2 DNA CONTROL	MBC024
AMPLIRUN® MYCOPLASMA GENITALIUM DNA CONTROL	MBC085
AMPLIRUN® MYCOPLASMA HOMINIS DNA CONTROL	MBC084
AMPLIRUN® NEISSERIA GONORRHOEAE DNA CONTROL	MBC075
AMPLIRUN® PAPILLOMAVIRUS TYPE 16 DNA CONTROL	MBC067
AMPLIRUN® PAPILLOMAVIRUS TYPE 18 DNA CONTROL	MBC068
AMPLIRUN® TREPONEMA PALLIDUM DNA CONTROL	MBC109
AMPLIRUN® TRICHOMONAS VAGINALIS DNA CONTROL	MBC079
AMPLIRUN® UREAPLASMA PARVUM DNA CONTROL	MBC133
AMPLIRUN® UREAPLASMA UREALYTICUM DNA CONTROL	MBC112
OTHER	
AMPLIRUN® CANDIDA ALBICANS DNA CONTROL	MBC010
AMPLIRUN® GARDNERELLA VAGINALIS DNA CONTROL	MBC108
AMPLIRUN® STREPTOCOCCUS AGALACTIAE DNA CONTROL	MBC071

AmpliRun® Total Extraction & Amplification Controls

Kit contents:

10 vials for single use that contain the complete microorganism in a matrix that mimics human specimens.

Description	Cat. No.
AMPLIRUN® TOTAL CT/NG CONTROL (EXUDATE) <i>Chlamydia trachomatis</i> & <i>Neisseria gonorrhoeae</i>	MBTC006
AMPLIRUN® TOTAL CT/NG CONTROL (URINE) <i>Chlamydia trachomatis</i> & <i>Neisseria gonorrhoeae</i>	MBTC003
AMPLIRUN® TOTAL CT/NG/TV/MGE CONTROL (SWAB) <i>Chlamydia trachomatis</i> , <i>Neisseria gonorrhoeae</i> , <i>Trichomonas vaginalis</i> , <i>Mycoplasma genitalium</i>	MBTC024
AMPLIRUN® TOTAL HSV-1/HSV-2/VZV CONTROL (EXUDATE) Herpes simplex virus & varicella-zoster virus	MBTC016